

---

# Project report

## SFC: Valves reconfiguration

*D Malan*

Department of Chemistry  
University of Pretoria  
6 December 2019

---

## Introduction

The aim of this project is reconfigure the valves of the SFC for resumption of chromatographic runs. The valves were un-configured for the column flow experiment. This was necessary because there is no CO<sub>2</sub> supply shut-off valve, so the system is always at CO<sub>2</sub> pressure. Because I wanted to disconnect and re-connect columns I needed to be able to shut down the CO<sub>2</sub> flow without depressurising the pump and losing all that expensive carbon dioxide.

---

## Experimental

The SFC pump was filled and set to a pressure of 200 atm. The flow through the restrictor was measured.

“The choice was a simple linear restrictor and we trusted that heating the end of the restrictor would prevent discrimination. The final linear restrictor was 800mm long and had an internal diameter of 0.050mm.”<sup>1</sup> “The SFC column we used in the SFC×GC system was ... HPLC columns (150mm × 4.6mm, 3µm particles) (Restek, Pinnacle DB Silica) connected in series.”<sup>2</sup>

The time taken for 50 ml of gaseous carbon dioxide to flow through a bubble flow meter was measured with a hand-held stopwatch, and typed into an Excel spreadsheet.

The flow was measured under a few configurations.

1. With the modifier valve reconfigured, in the LOAD position.
2. With the modifier valve in the RUN position.
3. With the modifier valve in the RUN position and the sample valve in the INJECT position.

---

## Observations

The flow is not as stable as one would wish. For the first two sets the flow was reliable and repeatable, but for the last set the flow started being a little bit puffy. Clipping the tip of the restrictor did not make a noticeable difference.

## Data analysis

Import the data:

```
In[1]:= data = Import[
  "\\chemserv\\CHEMBACK\\HPLC-GC\\Varian-SFC-x-GC#R2-26\\ANAL\\EgmontR\\Niel\\
  Dropbox\\2019_12_05_SFC_Flow.xlsx", "DataLegacy"];

In[11]:= data[[1]][[1 ;;]][[All, 3 ;;]] // TableForm

Out[11]//TableForm=
```

Column A	L1	Column B	L2	Column C	L3	Column D	L4	Column E
Fed via modifier valve and mixing chamber and sample valve.		Modifier valve in "Load"						
Volume of gas								
ml								
50.								
50.								
50.								
50.								
50.								
50.								
50.								

Column A	L1	Column B	L2	Column C	L3	Column D	L4	Column E
Fed via modifier valve and mixing chamber and sample valve.		Modifier valve in "Run"						
Volume of gas								
ml								
50.								
50.								
50.								
50.								
50.								
50.								
50.								
50.								

Column A	L1	Column B	L2	Column C	L3	Column D	L4	Column E
Fed via modifier valve and mixing chamber and sample valve.		Modifier valve in "Run"						
Volume of gas								
ml								
50.								
50.								
50.								
50.								
50.								
50.								
50.								
50.								
50.								

Copy & paste data from the cell above

```
In[14]:= data1 = { 429.79942693409737`
  452.4886877828054`
  438.59649122807014`
  436.04651162790697`
  443.1314623338257`
  436.04651162790697` } // Flatten

Out[14]= {429.799, 452.489, 438.596, 436.047, 443.131, 436.047}
```

In[15]:=

```
data2 = Flatten[
  {
    433.52601156069363`
    460.1226993865031`
    437.95620437956205`
    434.1534008683068`
    424.929178470255`
    408.7193460490463`
  }
]
```

Out[15]=

```
{433.526, 460.123, 437.956, 434.153, 424.929, 408.719}
```

In[42]:=

```
data3 = Flatten[
  {
    505.050505050505`
    406.5040650406504`
    407.6086956521739`
    404.3126684636119`
    402.1447721179625`
    410.39671682626545`
    406.5040650406504`
    412.08791208791206`
  }
]
```

Out[42]=

```
{505.051, 406.504, 407.609, 404.313, 402.145, 410.397, 406.504, 412.088}
```

In[29]:=

```
flow1 = Join[Split[data1], Table[1, Length[data1], 1], 2]
```

Out[29]=

```
{{429.799, 1}, {452.489, 1}, {438.596, 1}, {436.047, 1}, {443.131, 1}, {436.047, 1}}
```

In[31]:=

```
flow2 = Join[Split[data2], Table[2, Length[data2], 1], 2]
```

Out[31]=

```
{{433.526, 2}, {460.123, 2}, {437.956, 2}, {434.153, 2}, {424.929, 2}, {408.719, 2}}
```

In[43]:=

```
flow3 = Join[Split[data3], Table[3, Length[data3], 1], 2]
```

Out[43]=

```
{{505.051, 3}, {406.504, 3}, {407.609, 3}, {404.313, 3},
 {402.145, 3}, {410.397, 3}, {406.504, 3}, {412.088, 3}}
```

In[44]:=

```
flow = Join[flow1, flow2, flow3]
```

Out[44]=

```
{{429.799, 1}, {452.489, 1}, {438.596, 1}, {436.047, 1}, {443.131, 1},
 {436.047, 1}, {433.526, 2}, {460.123, 2}, {437.956, 2}, {434.153, 2},
 {424.929, 2}, {408.719, 2}, {505.051, 3}, {406.504, 3}, {407.609, 3},
 {404.313, 3}, {402.145, 3}, {410.397, 3}, {406.504, 3}, {412.088, 3}}
```

In[45]:=

```
flowavg = MapIndexed[ {Mean[#1], #2[[1]]} &, {data1, data2, data3}]
```


Out[45]=

```
{{439.352, 1}, {433.234, 2}, {419.326, 3}}
```

In[46]:=

```
ListPlot[{flow, flowavg},
  PlotStyle → {Blue, {Red, PointSize[0.02]}, ChartLabels → {0, 0, 1, 2}}]
```


Out[46]:=



In[47]:=


```
flowavgplot = ListPlot[flowavg, PlotStyle → {Red, PointSize[0.02]}]
```

Out[47]:=



```
In[51]:=
```

```
flowplot = ListPlot[flow]
```




Out[51]=

Combine the graphics:

```
In[54]:=
```

```
Show[flowplot
, flowavgplot

, PlotStyle → {Automatic, {Red, PointSize[0.09]}}
, Axes -> False
, Frame → {True, True, False, False}
, FrameLabel → {Style["Flow (ml/min)", Large, Black, Bold],
  Style["Run number", Large, Black, Bold]}
, FrameStyle → Directive[Medium, Bold, Black]
, ImageSize → 600
, PlotRange → {Automatic, Automatic}]
```



Out[54]=

---

## Results

The figure above summarizes the results. There is some variation in the flow for the different runs.

---

## Bibliography

- A. 1 Malan, Daniel. 2019. "Fast temperature programmed chromatography coupled to supercritical fluid chromatography (SFC×GC)."
- B. 2 Malan, Daniel. 2019. "Fast temperature programmed chromatography coupled to supercritical fluid chromatography (SFC×GC)."
- C. 3 <https://mathematica.stackexchange.com/questions/8997/cut-and-paste-data-from-a-spreadsheet>