Project report

SFC Flow

D Malan

Department of Chemistry
University of Pretoria
5 December 2019

Introduction

The aim of this project is to get a feeling for how much the flow through the SFC is determined by the number of columns.

Experimental

The SFC pump was filled and set to a pressure of 200 atm. The flow trough the restrictor was measured.

"The choice was a simple linear restrictor and we trusted that heating the end of the restrictor would prevent discrimination. The final linear restrictor was 800mm long and had an internal diameter of 0.050mm." 1

A different number of columns were connected in series, starting with one, and working up to 5. As a control, the flow without any columns connected were also measured.

"The SFC column we used in the SFC×GC system was ... HPLC columns (150mm × 4.6mm, 3μ m particles) (Restek, Pinnacle DB Silica) connected in series."²

The time taken for 50 ml of gaseous carbon dioxide to flow through a bubble flow meter was measured with a hand-held stopwatch, and typed into an Excel spreadsheet.


Data analysis


Import the data:

```
data = Import[
In[1]:=
           "\\\chemserv\\CHEMBACK\\HPLC-GC\\Varian-SFC-x-GC#R2-26\\ANAL\\EgmontR\\Nie1\\Dropbox
             \\2019_12_04_SFC_Flow.xlsx", "DataLegacy"];
        data[[2]][[1;; 10]] // Grid
In[ • ]:=
                   Time
                            Column A
                        Fed via Link 1
        4/12/2019
                         Volume of gas Time
                                                Time
                                                         Flow
                               m1
                                                min
                                                        ml/min ml/min 1.
                                         6.29 0.104833 476.948
                              50.
Out[ • ]=
                                         6.19 0.103167 484.653
                                         6.31 0.105167 475.436
                              50.
                              50.
                                         6.15 0.1025 487.805
                                         6.17 0.102833 486.224 482.213 1.
```

Data are in the following ranges. The list (a copy & paste from Excel³) gives the line number before the data starts, followed by the number of lines.

```
datarange = Flatten[#[[1]] + Range[#[[2]]] & /@
In[2]:=
             \{4, 5\}, \{13, 5\}, \{23, 5\}, \{34, 6\}, \{45, 6\}, \{56, 6\}, \{68, 6\}, \{79, 6\}, \{90, 6\}, \{102, 6\}\}\}
         {5, 6, 7, 8, 9, 14, 15, 16, 17, 18, 24, 25, 26, 27, 28, 35, 36, 37, 38,
Out[2]=
          39, 40, 46, 47, 48, 49, 50, 51, 57, 58, 59, 60, 61, 62, 69, 70, 71, 72, 73, 74,
          80, 81, 82, 83, 84, 85, 91, 92, 93, 94, 95, 96, 103, 104, 105, 106, 107, 108}
         Length[datarange]
In[ • ]:=
         57
Out[ • ]=
         flow = data[[2]][[datarange, {6, 8}]]
In[3]:=
         \{476.948, 1.\}, \{484.653, 1.\}, \{475.436, 1.\}, \{487.805, 1.\}, \{486.224, 1.\}, \{480., 2.\},
Out[3]=
          {483.871, 2.}, {483.871, 2.}, {487.013, 2.}, {477.707, 2.}, {473.186, 3.}, {465.839, 3.},
           {467.29, 3.}, {468.019, 3.}, {474.684, 3.}, {463.679, 4.}, {461.538, 4.}, {465.839, 4.},
          {468.75, 4.}, {475.436, 4.}, {464.396, 4.}, {455.235, 5.}, {465.116, 5.}, {460.829, 5.},
          \{473.934, 5.\}, \{472.441, 5.\}, \{464.396, 5.\}, \{449.775, 6.\}, \{459.418, 6.\},
          \{455.235, 6.\}, \{456.621, 6.\}, \{457.317, 6.\}, \{460.829, 6.\}, \{434.153, 7.\},
           \{436.681, 7.\}, \{432.9, 7.\}, \{436.681, 7.\}, \{436.047, 7.\}, \{430.416, 7.\}, \{465.116, 8.\},
          \{470.219, 8.\}, \{484.653, 8.\}, \{458.015, 8.\}, \{470.219, 8.\}, \{465.839, 8.\}, \{468.75, 9.\},
          {460.829, 9.}, {471.698, 9.}, {467.29, 9.}, {469.484, 9.}, {474.684, 9.}, {420.757, 10.},
          {424.328, 10.}, {421.348, 10.}, {415.512, 10.}, {421.941, 10.}, {424.929, 10.}}
       There are averages in the Excel sheet. The are at the last row of each range.
In[4]:=
         avgrange = Flatten[#[[1]] + #[[2]] & /@
             \{\{4, 5\}, \{13, 5\}, \{23, 5\}, \{34, 6\}, \{45, 6\}, \{56, 6\}, \{68, 6\}, \{79, 6\}, \{90, 6\}, \{102, 6\}\}\}
Out[4]=
         {9, 18, 28, 40, 51, 62, 74, 85, 96, 108}
In[5]:=
         flowavg = data[[2]][[avgrange, {7, 8}]]
         \{482.213, 1.\}, \{482.492, 2.\}, \{469.803, 3.\}, \{466.606, 4.\}, \{465.325, 5.\},
Out[5]=
          \{456.533, 6.\}, \{434.48, 7.\}, \{469.01, 8.\}, \{468.789, 9.\}, \{421.469, 10.\}\}
In[ • ]:=
         ListPlot[{flow, flowavg},
          PlotStyle \rightarrow {Blue, {Red, PointSize[0.02], ChartLabels \rightarrow {0, 0, 1, 2}}}
```





Prepare text for annotations:

```
textlist = \{"No. of columns", "1", "1", "2", "3", "3", "4", "5", "0", "0", "5"\}
In[8]:=
         {No. of columns, 1, 1, 2, 3, 3, 4, 5, 0, 0, 5}
Out[8]=
```

```
textgraphics =
In[12]:=
          MapIndexed[Graphics[Style[Text[#1, {495, #2[[1]] - 1}], Medium, Bold]] &, textlist];
```

Combine the graphics:

Show [flowplot In[14]:= , flowavgplot , textgraphics , PlotStyle → {Automatic, {Red, PointSize[0.09]}} , Axes -> False , Frame → {True, True, False, False} , FrameLabel → {Style["Flow (ml/min)", Large, Black, Bold], Style["Run number", Large, Black, Bold]} , FrameStyle → Directive [Medium, Bold, Black] , ImageSize → 600 , PlotRange → { {414, 500}, Automatic}] 5 Out[14]= No. of columns 440 460 480 420 Flow (ml/min)

Results

The figure above summarizes the results. It is clear more columns leads to lower flow, but that inter-run repeatability is not very high. For runs 6 and 10 all five columns were in series, but there is a significant difference in flow. For runs 1 and two, with only one column connected, the flow is higher than the flow for runs 8 and 9, where there were no columns. None of these observations suggest that the SFC separations will be significantly impacted with 5 columns in series.

Bibliography

- A. 1 Malan, Daniel. 2019. "Fast temperature programmed chromatography coupled to supercritical fluid chromatography (SFC×GC)."
- B. 2 Malan, Daniel. 2019. "Fast temperature programmed chromatography coupled to supercritica fluid chromatography (SFC×GC)."
- C. 3 https://mathematica.stackexchange.com/questions/8997/cut-and-paste-data-from-a-spreadsheet